





















NYLON - Monofilament Polyamide Suture

NEEDLE			MATERIAL				
Length and description			Length	Gauge			
				6/0	5/0	4/0	3/0
				0.7	1	1.5	2
15mm							
Curved cutting fine	▲		75cm			LM3251	
Curved cutting fine	▲		50cm	LM3249			
16mm							
Curved cutting	▲		35cm	LM3359			
Curved cutting	▲		50cm	LM3369	LM3360	LM3361	LM3362
Curved reverse cutting	▼		45cm				LM3362R
19mm							
Curved reverse cutting	▼		45cm		LM3260R	LM3261R	LM3262R
Curved cutting	▲		45cm		LM3330	LM3331	LM3332
20mm							
Curved cutting fine	▲		45cm			LM3266	
25mm							
Curved cutting fine	▲		75cm			LM3269	








NYLON - Monofilament Polyamide Suture

NEEDLE			MATERIAL				
Length and description			Length	Gauge			
				2/0	0	1	2
				3	3.5	4	5
26mm							
Curved cutting	▲		50cm	LM3282			
Curved reverse cutting	▼		45cm	LM3283R			
30mm							
Half circle round bodied	●		1.5m LOOP	LM3013L	LM3014L	LM3015L	
Curved cutting	▲		75cm	LM3343	LM3344		
Curved reverse cutting	▼		50cm	LM3003R			
40mm							
Half circle round bodied	●		1m	LM3033	LM3034	LM3035	LM3036
Half circle round bodied	●		1.5m LOOP	LM3033L	LM3034L	LM3035L	LM3036L
Curved cutting	▲		45cm	LM3403	LM3404		
45mm							
Curved reverse cutting	▼		1m	LM3496R	LM3497R		



NYLON - Monofilament Polyamide Suture

NEEDLE		MATERIAL					
Length and description	Length	Gauge					
		2/0	0	1	2		
		3	3.5	4	5		
50mm							
Half circle round bodied	●		1 m		LM3064		
Half circle round bodied	●		1.5 m LOOP		LM3064L		
Half circle round bodied heavy	●		1 m			LM3065H	LM3066H
60mm							
Straight cutting	▲		75cm	LM3453	LM3454		
90mm							
Curved cutting	▲		1 m	LM3483	LM3484	LM3485	LM3486



NYLON - Monofilament Polyamide Suture

TECHNICAL INFORMATION

STERILE MONOFILAMENT NON - ABSORBABLE SUTURE OF POLYAMIDE - BLUE

- ❖ Strong uniform monofilament strand
- ❖ Smooth surface does not support bacterial growth
- ❖ Unformed elasticity permits its use in ultra fins sizes in all types of microsurgeries

Nylon sutures are a monofilament, synthetic, non-absorbable sterile surgical suture composed of polyamide 6, 6/6. Nylon sutures are dyed blue with pigment C174160 pigment blue, 1.5 irgalite Blue BLP- Phthalocyaninato (2) Copper.

Nylon sutures complies with the requirements of the European Pharmacopoeia for sterile Polyamide 6, 6/6 sutures, the United States Pharmacopoeia for non-absorbable surgical sutures and with the essential requirements of the Medical Device Directive 93/42/EEC.

CE0088 applies to our Nylon sutures.

Indications

Nylon sutures are intended for use in general skin and fascia closure procedures.

Application

Sutures should be selected and implanted depending on surgical technique, surgical experience, patient's condition and wound size.

Sterilised

By either Gamma Irradiation or Ethylene Oxide