












SILK – Braided Silk Suture

NEEDLE			MATERIAL				
Length and description			Length	Gauge			
				6/0	5/0	4/0	3/0
				0.7	1	1.5	2
13mm							
Curved round bodied	●	⤿	75cm	LM2559			
15mm							
Curved cutting fine	▲	⤿	45cm	LM2769	LM2770	LM2771	LM2772
16mm							
Curved round bodied	●	⤿	75cm		LM2790	LM2791	
Half circle round bodied	●	⤿	75cm		LM2800	LM2801	LM2802
19mm							
Curved reverse cutting	▼	⤿	45cm		LM2900R	LM2901R	
Curved cutting	▲	⤿	45cm		LM2900		
20mm							
Half circle round body	●	⤿	75cm		LM2920	LM2921	LM2922
Curved round bodied	●	⤿	75cm				LM2912
Curved cutting fine	▲	⤿	45cm		LM2905	LM2906	








SILK – Braided Silk Suture

NEEDLE			MATERIAL				
Length and description	Length	Gauge					
		3/0	2/0	0	1		
		2	3	3.5	4		
22mm							
Half circle cutting	▲		45cm	ST 577	ST 578		
25mm							
Half circle round bodied	●		75cm	LM2942	LM2943	LM2944	
Curved cutting fine	▲		45cm	LM2972	LM2973		
30mm							
Half circle round bodied	●		75cm	LM2952	LM2953	LM2954	LM2955
Curved cutting	▲		75cm	LM2812	LM2813	LM2814	LM2815
35mm							
Half circle round bodied	●		75cm		LM2993	LM2994	LM2995
Curved reverse cutting	▼		75cm	LM2562R	LM2563R	LM2564R	LM2565R
Curved cutting	▲		75cm	LM2562			
40mm							
Curved cutting heavy	▲		75cm		LM2883		



SILK – Braided Silk Suture

NEEDLE			MATERIAL				
Length and description	Length	Gauge					
		2/0	0	1	2		
		3	3.5	4	5		
40mm							
Half circle round bodied	●		75cm	LM2513	LM2514	LM2515	LM2516
Half circle round bodied heavy	●		75cm	LM2513H	LM2514H	LM2515H	LM2516H
Half circle reverse cutting	▼		75cm		LM2524R	LM2525R	LM2526R
45mm							
Curved reverse cutting	▼		75cm	LM2583R	LM2584R		
55mm							
Straight cutting	▲		1m	LM2533	LM2534		



SILK – Braided Silk Suture

NON – NEEDLE		LENGTH		GAUGE			
		4/0	3/0	2/0	0	1	2
		1.5	2	3	3.5	4	5
Cut lengths	13 x 60cm	LM2441	LM2442	LM2443	LM2444	LM2445	LM2446
Cut lengths	1.8m	LM2501	LM2502	LM2503	LM2504	LM2505	LM2506
Non – Sterile - Reels	22m	LR2491	LR2492	LR2493	LR2494	LR2495	LR2496
Non – Sterile - Reels	90m	LR2601	LR2602	LR2603	LR2604	LR2605	LR2606



SILK – Braided Silk Suture

TECHNICAL INFORMATION

STERILE BRAIDED NON - ABSORBABLE SUTURE OF NON – CAPILLARY SILK – BLACK AND WHITE

- ❖ 100% natural braided silk
- ❖ Developed by a special purification process (degumming). This process ensures a compact braid while allowing the filaments to retain their natural body and elasticity
- ❖ Specially treated to make the thread serum proof and water repellent
- ❖ Does not swell in tissue fluids
- ❖ Easy to handle and ensures perfect knot security
- ❖ Better visibility due to the black colour

Silk suture are a braided, non-absorbable sterile surgical suture composed of silk fibrils. This is derived from the domestic species *Bombyx Mori* (B> mori) of the family Bombycidae. Silk for braided material is impregnated with either wax or silicone and is dyed black with C.153185 Sulphol Black 1.

Silk sutures complies with the requirements of the European Pharmacopoeia for silk sutures, the United States Pharmacopoeia for non - absorbable surgical sutures and with the essential requirements of the Medical Device Directive 93/42/EEC.

CE0088 applies to our Silk sutures.

Indications

Silk sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures

Application

Sutures should be selected and implanted depending on surgical technique, surgical experience, patient's condition and wound size.

Sterilised

By either Gamma Irradiation or Ethylene Oxide